



WASP[®] IN THE WORLD THE LARGEST STREAKER INSTALLED BASE

Please contact your local bioMérieux representative for availability of these solutions in your country.

Manufactured by Via F. Perotti 10 • 25125 Brescia • Italy



Distributed by bioMérieux S.A. • 69280 Marcy l'Etoile • France • Tel.: + 33 (0)4 78 87 20 00 • Fax: +33 (0)4 78 87 20 90 www.biomerieux.com

WASPLab[®], AVAILABLE FOR YOUR LAB TOMORROW



Your need for extra productivity may evolve over time. No problem! Count on us to find the right set-up to suit your new needs. Your WASP[®] can be a launch pad for increased lab efficiency by evolving to WASPLab[®] over time.

It's a smart investment in the age of the automated lab.



WASP[®] EMPOWERED BY BIOMÉRIEUX We adapt to you. You don't have to adapt to us.





THE AGILE & FLEXIBLE AUTOMATED SPECIMEN PROCESSING SYSTEM



YOUR LAB. YOUR WAY. YOU

Build your workstation to suit your lab! WASP[®] is a **unique** needs change. Combine it with services available only through **V**

THE MOST AUTON

WASP° fits into your lab

Automated Gram slide preparation



Pipetor module from 10 µl to 900 µl



Automated broth preparation



QUALITY FROM START TO FINISH

As microbiologists, we understand the need to pamper specimens to get the best results. WASP® offers easy solutions for quality colony growth:

- Use with compatible Liquid Based Media (LBM[°]) and Liquid Amies
 Collection SWAB (ESWAB[°]) ranges for specimen protection
 throughout transport, optimized sample conservation & release and
 streamlined workflow with standardized sample containers
- Fully validated compatibility with **bioMérieux's broad range of high** quality pre-plated media (PPM)

CUSTOMIZATION FOR OPTIMAL COLONY ISOLATION

Keep the standardization of automation and team it with specimen inoculation flexibility to suit your lab needs.

- Total versatility with an extensive and customizable library of classic streaking patterns that reproduce human gestures
- Choose your specimen transfer volumes: Loops for 1 μl, 10 μl or 30 μl and even twin loops for efficient biplate streaking
- Optional **loops sterilization between quadrants** for optimal colony isolation of samples with high bacterial loads
- Work with **any specimen type** (no pipetting required)

R WASP[®].

open-platform, modular system that grows with you as your NASP[®] empowered by bioMerieux for ultimate automation.

IATED SPECIMEN PROCESSOR AVAILABLE

and lets you pick and mix the automation add-ons you need – we adapt to you, you don't have to adapt to us!



Continuous Loading



Urine Turntable



Antibiotic Disc Deposit system

Why choose WASP[®] empowered by bioMérieux ?

WASP[®] empowered by bioMérieux is, simply, more. We back your WASP[®] instrument with **flexible**, **innovative services** developed using our 50 years of microbiology expertise.

CONCRETELY, WE:

- Partner with you for the life of your WASP through our **5-star support service package**
- Provide relevant training and education by microbiology experts
- Leverage our **microbiology know-how** to optimize your WASP and its use
- Smoothly integrate automation and informatics facilitating the continuous improvement of your processes

